

# Healthy Gas Permeable Contact Lens Habits


Copyright © 2012 AOCLE. All Rights Reserved.

## Congratulations on your new GP contact lenses!

### To ensure continued success with your GP lenses, review these healthy lens habits.

1. Wash your hands thoroughly with soap, rinse and dry them before handling your lenses.<sup>1,2</sup>



2.  Do not use tap water to clean or handle your lenses or to clean the case.<sup>3-7</sup>


3. Your solution has been chosen specifically for your type of lens. Do not change without discussing this with your eye doctor.<sup>8-11</sup>  
In addition, contact lens solutions should not be used after the expiration date.<sup>12</sup>  
Your chosen solution is \_\_\_\_\_.

4. To maintain comfortable lens wear, rub your lenses with the prescribed solution to remove protein, oil and make-up.<sup>13-17</sup>

5. Store brand solutions are often old formulations of solutions and may not be compatible with your type of lens.<sup>18-21</sup>

6. Always recap your solution bottle.<sup>1,22</sup>



7.  Do you want to nap or sleep with your lenses? Ask your doctor if your lenses are designed for that.<sup>23,24</sup> Sleeping in your contact lenses, if not prescribed for overnight wear, increases your risk for eye infection.<sup>25-29</sup>

8. You should store your spare pair dry in the case until needed. Then clean the lens(es) and store in disinfecting solution overnight prior to application.<sup>30</sup>

9. Before using any type of eye drop (medication or artificial tear), ask your doctor if the lenses need to be removed.<sup>31-33</sup>



10. Clean the case with your prescribed solution and let air dry completely, uncapped and upside down, on a paper towel. It is recommended to replace your case at least every 3 months.<sup>2, 34-38</sup>



11. Empty completely the solution in the lens case every day. Never add more solution if there is solution already there.<sup>16, 35, 39, 40</sup>



12. Healthy lens wear should not induce redness, discomfort, or visual disturbance. If you have any concerns remove your lenses and consult your eye care professional.<sup>41, 42</sup>

13. Your eye doctor has selected a wearing and replacement schedule for you. These are important to be followed.<sup>43-45</sup>

Your next contact lens checkup is:

\_\_\_\_\_

Your contact lenses must be replaced every:

3 months  9 months  1 year  as needed

# References

1. Szczotka-Flynn LB, Pearlman E, et al. Microbial contamination of contact lenses, lens care solutions, and their accessories: a literature review. *Eye & Contact Lens* 2010;36(2):116-129.
2. Wu Y, Carnt N, et al. Contact lens user profile, attitudes and level of compliance to lens care. *Cont Lens Ant Eye* 2010;33(4):183-188.
3. Shoff ME, Joslin CE, et al. Efficacy of contact lens systems against recent clinical and tap water acanthamoeba isolates. *Cornea* 2008;27(6):713-719.
4. Joslin CE, Tu EY, et al. The association of contact lens solution use and Acanthamoeba keratitis. *Am J Ophthalmol* 2007;144(2):169-180.
5. Booton G, Joslin CE, et al. Genotypic identification of acanthamoeba sp. Isolates associated with an outbreak of acanthamoeba keratitis. *Cornea* 2009;28(6):673-676.
6. Amos CF, George MD. Clinical and laboratory testing of silver-impregnated lens case. *Cont Lens Ant Eye* 2006;29(5):247-255.
7. Hall BJ, Jones L. Contact lens cases: the missing link in contact lens safety? *Eye & Contact Lens* 2010;36(2):101-105.
8. Keir A, Woods CA, et al. Clinical performance of different care systems with silicone hydrogel contact lenses. *Cont Lens Ant Eye* 2010;33:189-195.
9. Carnt N, Evans V, et al. IER Matrix update: adding another silicone hydrogel. *Contact Lens Spectrum* 2008;23(3):4-43.
10. Andrasko G, Ryan K. Corneal staining and comfort observed with traditional and silicone hydrogel lenses and multipurpose solution combinations. *Optometry* 2008;79(8):444-454.
11. Jones L, McDougail N, et al. Asymptomatic corneal staining and comfort observed with the use of balafilcon silicone-hydrogel contact lenses disinfected with a polyaminopropyl biguanide-preserved care regimen. *Optom Vis Sci* 2002;79(12):753-761.
12. [www.fda.gov/NewsEvents/Newsroom/PressAnnouncements/2009/ucm167187.htm](http://www.fda.gov/NewsEvents/Newsroom/PressAnnouncements/2009/ucm167187.htm).
13. Sluh KL, Hu J, Sibley M. The microbiological benefit of cleaning and rinsing contact lenses. *Int Contact Lens Clin* 1985;12:135-142 (check pg numbers).
14. Henry VA. The role of compliance. *Rev Cornea & Contact Lenses*, 2009(10):27-32.
15. Rosenthal RA, Henry CL, Schlech BA. Contribution of regimen steps to disinfection of hydrophilic contact lenses. *Contact Lens & Anterior Eye* 2004;27(3):149-156.
16. Butcko V, McMahon TT, Joslin CE, Jones L. Microbial keratitis and the role of rub and rinsing. *Eye & Contact Lens* 2007;33(6):421-423.
17. Cho P, Cheng SY, Chan WY, Yip WK. Soft contact lens cleaning: rub or no rub? *Ophthalmol Physiol Opt* 2009;29(1):49-57.
18. Wu Y, Carnt N, et al. Contact lens and lens storage case cleaning instructions: Whose advice should we follow? *Eye & Contact Lens* 2010;36(2):68-72.
19. Gromacki S. The truth about generics. *Contact Lens Spectrum* 2005;20(12):24.
20. Gromacki S. Solution confusion. *Rev of Cornea & Contact Lenses (Rev Optom suppl)*;147(1).
21. Ward MA. How private label solutions affect your practice. *Contact Lens Spectrum* 2006;21(3):25.
22. Collins M, Coulson J, et al. Contamination of disinfection solution bottles used by contact lens wearers. *CLAO J* 1994;20(1):32-36.
23. Cho P, Cheung S. Napping in contact lenses-an underestimated risk. *Contact Lens Spectrum* 2010;25(2).
24. Dumbleton K. Adverse events with silicone hydrogel continuous wear. *Cont Lens Ant Eye* 2002;25:137-146.
25. Schein, O.D., et al., The relative risk of ulcerative keratitis among users of daily-wear and extended-wear soft contact lenses. A case-control study. *Microbial Keratitis Study Group. N Engl J Med* 1989;321(12):773-8.
26. Poggio, E.C., et al., The incidence of ulcerative keratitis among users of daily-wear and extended-wear soft contact lenses. *N Engl J Med* 1989;321(12): 779-83.
27. Cheng, K.H., et al., Incidence of contact-lens-associated microbial keratitis and its related morbidity. *Lancet* 1999;354(9174):181-5.
28. Schein, O.D., et al., The incidence of microbial keratitis among wearers of a 30-day silicone hydrogel extended-wear contact lens. *Ophthalmology* 2005; 112(12):2172-9.
29. Dart JK, Radford CF, Minassian D, Verma S, Stapleton F. Risk factors for microbial keratitis with contemporary contact lenses: a case-control study. *Ophthalmology* 2008;115:1647-54.
30. Bennett ES, Wagner H. Rigid lens care and patient education. In Bennett ES, Weissman BA. *Clinical Contact Lens Practice*. Philadelphia, Lippincott Williams & Wilkins, 2005:277-294.
31. Tabuchi N, Watanabe T, et al. Adsorption of actives in ophthalmological drugs for over-the-counter on soft contact lens surfaces. *J Oleo Sci* 2009;58(1):43-52.
32. Kading D. A two-week clinical evaluation of the safety of systane ultra in contact lens-wearing patients. *Clin Ophthalmol* 2010;4:27-32.
33. Nichols J, Sinnott L, et al. Hydrogel contact lens binding induced by contact lens rewetting drops. *Optom Vis Sci* 2008;85(4):236-240.
34. Yung A, Boost M, et al. The effect of compliance enhancement strategy (self-review) on the level of lens care compliance and contamination of contact lenses and lens care accessories. *Clin Exp Optom* 2007;90(3):190-202.
35. [www.fda.gov/MedicalDevices/ProductsandMedicalProcedure/HomeHealthandConsumer/ConsumerProducts/Contactlenses/](http://www.fda.gov/MedicalDevices/ProductsandMedicalProcedure/HomeHealthandConsumer/ConsumerProducts/Contactlenses/)
36. Willcox M, Carnt N, et al. Contact lens case contamination during daily wear of silicone hydrogels. *Optom Vis Sci* 2010;87:456-464.
37. Wu Y, Zhu H, et al. Impact of air-drying lens cases in various locations and positions. *Optom Vis Sci* 2010;87:465-468.
38. Wu YT, Zhu H, Willcox M, Stapleton F. Removal of biofilm from contact lens storage cases. *Invest Ophthalmol Vis Sci* 2010;51(12):6329-6333.
39. Proenca-Pina J, Ssi Yan Kai I, et al. Fusarium keratitis and endophthalmitis associated with contact lens wear. *Int Ophthalmol* 2010;30(1):131-139.
40. Kilvington S, Powell CH, Lam A, Lonnen J. Antimicrobial efficacy of multi-purpose contact lens disinfectant solutions following evaporation. *Cont Lens Ant Eye* 2011;Mar 8 (Epub ahead of print).
41. FDA Consumer Health Information. "Focusing on contact lens safety" [www.fda.gov/consumers/features/contactlens1107.html](http://www.fda.gov/consumers/features/contactlens1107.html). Accessed 4/08/2011.
42. Kame RT. Clinical management of hydrogel-induced edema. *Am J Optom Physiol Opt* 1976;15:61-63.
43. Dumbleton K, Richter D, et al. Compliance with contact lens replacement in Canada and the United States. *Optom Vis Sci* 2010;87(2):131-139.
44. Dumbleton K, Woods C, et al. Comfort and vision with silicone hydrogel lenses: effect of compliance. *Optom Vis Sci* 2010;87(6): 421-425.
45. US Food and Drug Administration: [www.fda.gov/medicaldevices/productsandmedicalprocedures/homehealthandconsumer/consumerproducts/contactlenses/default.htm](http://www.fda.gov/medicaldevices/productsandmedicalprocedures/homehealthandconsumer/consumerproducts/contactlenses/default.htm). Accessed 2/20/2011.